



**MILA**   
INTERNATIONAL, INC.  
*Medical Instrumentation for Animals*

**2012 Catalog**




**MILA International** is an innovative company which has developed veterinary medical products otherwise commercially unavailable. We appreciate the opportunity to serve the world's diverse veterinary community.

Accessories & Components	34-36	Fluid Control Dial	30
Administration Sets	30-31	Foley Catheters	15
Airway Management	25	Gastrostomy Devices	18
Antimicrobial Catheters	4-7	GDV / Bloat Tube	23
Arthroscopy Fluid Line	31	Grenade - Suction Bulbs	20 & 22
Barbed Adaptors	34	Guidewire Catheters	6
Blake Drain	22	Guidewires	8 & 34
Bloat Tube	23	Guttural Pouch	28
Blood Collection	30	Heimlich Valve	21
Bone Marrow Aspiration Needle	19	Infusion Pumps	
Broncho Alveolar Lavage - Equine	32	Clear Cuff Pressure Infusors	31
Broncho Alveolar Lavage - Small Animal	23	Elastomeric Administration Pump	26
Caps	13	Infu-Disk	27
Catheters		Spring Powered Syringe Pump	26
BAL	23/32	Insemination Catheter	28
Diffusion / Wound	24	Instructional DVD	13
Endoscopic Aspiration & Delivery Catheters	23/28	IV Backpack	32
Epidural Catheters	24	Jackson Pratt Drains	22
Foley	15	Jejunostomy Tubes	18
Guttural Pouch	28	Joint Infusion System	33
Pain Management	24	Lance-A-Lot	28
Tomcat	14	Long Term Catheters	4-8
Urinary Catheters	14-15	MicroClave	12
Venous Access		Multi-purpose Drain	22
Antimicrobial Long Term Catheter	4-7	Nasogastric Feeding Tubes	16
Equine CVP	10	Nasolacrimal Duct Catheter	27
Extended Use MILACATH	9	Needles	19
Guidewire Introducer Technique	5	Oocyte Collection Set	29
Long Line Kit	10	Oxygen Tubes	16
PICC / Silicone Catheter	11	Pain Management	
Peel Away Introducer Technique	5	Diffusion / Wound Catheter	24
Polyurethane - Long Term Catheter	6-7	Epidural Kit	24
Short Term	9-10	Peel Away Introducers	8
Introcan Safety	10	PEG Kit	18
VeinCath	11	Pinch Clamps	36
Wound / Pain Management	24	Plasmapheresis Catheter	31
Centesis Valve	21	Pressure Infusors	31
Chest Tubes	20	Reproduction	29
Christmas Trees	34	Short Term Catheters	9-10
Connectors	34-36	Spinal Needles	19
Controlled Flushing Set	19	SPL	27
Conversion Chart	36	Stainless Steel Introducer	32
Cubby - Low Profile Gastrostomy Device	17	Stomach Tube	23
CVP Manometer	11	Stopcock	35
Dilator	8 & 34	Subpalpebral Eye Lavage	27
Diffusion Catheter	24	Suction Connector	21
Distributors	38	Surcingle	32
Drains	22	Suture Wings	36
Drain Adaptors	21-22	Technical Notes - Catheter Insertion Methods	5
DVD	13	Tendon Sheath	33
Embryo Transfer	29	Thoracic Catheters	20
Endoscope Products	23 & 28	Tomcat Urinary Catheter	15
Endotracheal Tubes	25	Trans-Tracheal Wash Kit	23 & 32
Enteral Feeding Products	16-19	TTL	
Epidural Kit	24	Small Animal	23
Esophagostomy Feeding Tubes	17	Large Animal	32
Extended Use MILACATH	9	Tuohy Needles	19
Extension Tubes	12	Urinary Catheters	15-16
Eye Lavage Kit	27	Uterine Flush Kit	29
Ferret Urinary Catheter	15	Veincath	11
Filter	35	Wound Catheter	24
Fluid Collection & Delivery	30-31	Wound Flushing Line	31


## Long Term Polyurethane Catheter Kits

- Non-reactive Polyurethane.
- Available with guidewire introducer or peel away introducer.
- Single, double, or triple lumen catheters for long term use. Can be left indwelling 3-4 weeks.
- Use for pressure monitoring, blood sampling, TPN and administration of drugs or fluids.
- Normally closed valve helps prevent loss of catheter due to accidental disconnect.
- Antimicrobial catheters are coated on both internal and external surfaces with silver sulfadiazine for the immune suppressed patient. Silver Sulfadiazine coating provides a precautionary protective enhancement.


**Single Lumen Catheters** - Available with integrated or detached extension tube. The 14Ga catheter has a side and distal exit port.




**Multi-Lumen Catheters** - Separate lumens with separate exit points.



Double Lumen



Triple Lumen




Quad Lumen

Long term catheter insertion methods include guidewire introduction or peel away introduction. Procedures for each method and catheter sizes are detailed on the following pages.



Long Term MILACATHs



Guidewire Introduction



Peel Away Introduction

## GUIDEWIRE INTRODUCTION of a Long Term MILACATH®



1 Photo shows introducer being inserted into vein.



2 Guidewire being advanced from housing. Once guidewire advanced to desired length, remove housing and introducer needle while holding onto the guidewire. **ALWAYS HOLD ONTO GUIDEWIRE.**



3 Slide the distal end of the catheter onto guidewire and advance until guidewire comes through the catheter hub. Grasp wire above the catheter hub. Advance catheter into vein by holding catheter near insertion site and pushing into place. Remove guidewire.



4 Suture catheter individually at each wing toward insertion site and across suture groove.



↑ GUIDEWIRE INTRODUCER KIT includes choice of introducers, guidewire, catheter, suture wings, MicroClave and optional dilator.

→ PEEL AWAY INTRODUCER KIT includes peel away catheter-over-needle introducer, catheter, suture wings, normally closed valve and tethered cap.

**For insertion demonstration, order DVD100**

## PEEL AWAY INTRODUCTION of a Long Term MILACATH®



1 Hold peel away introducer, bevel up, directly above vein with index finger firmly placed on top of stylet. Advance introducer into vessel far enough to be sure it is well seated.



2 Good flashback should be noted or blood easily drawn into syringe through stylet. Remove stylet while stabilizing outer portion of introducer.



3 Advance long term catheter through catheter of peel away introducer.



4 Hold catheter hub or extension set against animal. Grasp blue tabs with thumb and fore-finger of both hands. Gently pull up and out to separate blue tabs.



5 Introducer should split above the skin.



6 The catheter will back out of insertion site an inch or so during peel away procedure but can be easily pushed back.



7 Secure catheter individually at each wing toward insertion site and at suture groove. Attach additional suture wing to extension tube and suture to secure.

*Photos courtesy of Hagyard Equine Medical Institute*

## Long Term MILACATH Kits using GUIDEWIRE INTRODUCTION

Single Lumen			
DESCRIPTION		ITEM #	PRICE
19Ga (3.5Fr) x 15cm (6in)	20Ga catheter introducer, dilator	1906	\$20.85
19Ga (3.5Fr) x 25cm (10in)	20Ga catheter introducer, dilator	1910	\$20.85
16Ga (5Fr) x 15cm (6in)	18Ga catheter introducer, dilator	SA1606	\$20.85
16Ga (5Fr) x 15cm (6in)	with integrated extension set, choice of 16Ga & 18Ga introducers, dilator	1606	\$26.95
16Ga (5Fr) x 15cm (6in)	extension set, choice of 16Ga & 18Ga introducers, dilator	1606-2P	\$24.20
16Ga (5Fr) x 25cm (10in)	18Ga catheter introducer, dilator	SA1610	\$20.85
16Ga (5Fr) x 20cm (8in)	with integrated extension set, choice of 16Ga & 18Ga introducers, dilator	1610	\$26.95
16Ga (5Fr) x 20cm (8in)	antimicrobial catheter, choice of introducers, dilator	A1610	\$38.60
16Ga (5Fr) x 20cm (8in)	extension set, choice of 16Ga & 18Ga introducers, dilator	1610-2P	\$24.20
16Ga (5Fr) x 25cm (10in)	with integrated extension set, choice of 16Ga & 18Ga introducers, dilator	1625	\$26.95
16Ga (5Fr) x 25cm (10in)	extension set, choice of 16Ga & 18Ga introducers, dilator	1625-2P	\$24.20
14Ga (6Fr) x 15cm (6in)	with integrated extension set, choice of 14Ga & 16Ga introducers, dilator	1406	\$26.95
14Ga (6Fr) x 15cm (6in)	extension set, choice of 14Ga & 16Ga introducers, dilator	1406-2P	\$24.20
14Ga (6Fr) x 20cm (8in)	with integrated extension set, choice of 14Ga & 16Ga introducers, dilator	1410	\$26.95
14Ga (6Fr) x 20cm (8in)	antimicrobial catheter, choice of introducers, dilator	A1410	\$38.60
14Ga (6Fr) x 20cm (8in)	extension set, choice of 14Ga & 16Ga introducers, dilator	1410-2P	\$24.20
12Ga (8Fr) x 20cm (8in)	with integrated extension set, choice of 14Ga introducers, dilator	1210	\$28.00

Double Lumen			
4Fr x 8cm (3.25in)	with two 22Ga lumens, 18Ga catheter introducer, dilator	1540	\$43.25
4Fr x 13cm (5.25in)	with two 22Ga lumens, 18Ga catheter introducer, dilator	1500	\$43.25
5Fr x 8cm (3.25in)	with a 18Ga & 20Ga lumen, 16Ga catheter introducer, dilator	1640	\$43.25
5Fr x 13cm (5.25in)	with a 18Ga & 20Ga lumen, 16Ga catheter introducer, dilator	1600	\$43.25
5Fr x 25cm (10in)	with a 18Ga & 20Ga lumen, 16Ga catheter introducer, dilator	1650	\$43.25
7Fr x 20cm (8in)	with 14Ga & 18Ga lumen, 14Ga & 16Ga introducers, dilator	1620	\$43.25
7Fr x 20cm (8in)	antimicrobial catheter, choice of introducers, dilator	A1620	\$55.10
7Fr x 30cm (12in)	with a 14Ga & 18Ga lumen, choice of introducers, dilator	1630	\$43.25
7Fr x 60cm (24in)	with two 16Ga lumens, 14Ga & 16Ga introducers, dilator	1660	\$43.25

Triple Lumen			
5.5Fr x 8cm (3.25in)	with one 20Ga & two 22Ga lumens, 16Ga catheter introducer, dilator	1740	\$46.30
5.5Fr x 13cm (5.25in)	with one 20Ga & two 22Ga lumens, 16Ga catheter introducer, dilator	1700	\$46.30
7Fr x 15cm (6in)	with one 16Ga & two 18Ga lumens, 14Ga & 16Ga introducers, dilator	1720	\$46.30
7Fr x 60cm (24in)	with one 16Ga & two 18Ga lumens, 14Ga & 16Ga introducers, dilator	1760	\$46.30

Quadruple Lumen			
8.5Fr x 20cm (8in)	with one 16Ga, one 14Ga and two 18Ga lumens, introducers, dilator	1820	\$ 48.50

Slide clamps (#8599) sold separately, see page 12.  
Guidewires sold separately, see page 8.

GUIDEWIRE INTRODUCER  
KIT includes choice of introducers,  
guidewire, catheter, suture wings,  
MicroClave and optional dilator.



## Long Term MILACATH Kits using PEEL AWAY INTRODUCTION

### Single Lumen

DESCRIPTION		ITEM #	PRICE
19Ga (3.5Fr) x 15cm (6in)	4Fr x 1.5 inch peel away introducer	PI-1906	\$20.85
19Ga (3.5Fr) x 25cm (10in)	4Fr x 1.5 inch peel away introducer	PI-1910	\$20.85
16Ga (5Fr) x 15cm (6in)	5Fr x 2.5 inch peel away introducer	PI-506	\$20.85
16Ga (5Fr) x 25cm (10in)	5Fr x 2.5 inch peel away introducer	PI-510	\$20.85
16Ga (5Fr) x 15cm (6in)	with integrated extension set, 5Fr x 1.5 inch peel away introducer	PI-1606	\$26.95
16Ga (5Fr) x 15cm (6in)	extension set, 5Fr x 1.5 inch peel away introducer	PI-1606-2P	\$24.20
16Ga (5Fr) x 20cm (8in)	with integrated extension set, 5Fr x 2.5 inch peel away introducer	PI-1610	\$26.95
16Ga (5Fr) x 20cm (8in)	antimicrobial with integrated extension set, 5Fr x 1.5 inch peel away introducer	API-1610	\$38.60
16Ga (5Fr) x 20cm (8in)	extension set, 5Fr x 1.5 inch peel away introducer	PI-1610-2P	\$24.20
16Ga (5Fr) x 60cm (24in)	5Fr x 2.5 inch peel away introducer	PI-1624	\$20.85
16Ga (5Fr) x 25cm (10in)	with integrated extension set, 5Fr x 2.5 inch peel away introducer	PI-1625	\$25.95
16Ga (5Fr) x 25cm (10in)	extension set, 5Fr x 2.5 inch peel away introducer	PI-1625-2P	\$24.20
16Ga (5Fr) x 120cm (48in)	5Fr x 2.5 inch peel away introducer	PI-1648	\$20.85
14Ga (6Fr) x 20cm (8in)	with integrated extension set, 6Fr x 1.5 inch peel away introducer	PI-1410	\$26.95
14Ga (6Fr) x 20cm (8in)	antimicrobial with integrated extension set, 6Fr x 1.5 inch peel away introducer	API-1410	\$38.60
14Ga (6Fr) x 20cm (8in)	extension set, 6Fr x 1.5 inch peel away introducer	PI-1410-2P	\$24.20

### Double Lumen

4Fr x 8cm (3.25in)	with two 22Ga lumens, 4Fr x 1.5 inch peel away introducer	PI-1540	\$43.25
4Fr x 13cm (5.25in)	with two 22Ga lumens, 4Fr x 1.5 inch peel away introducer	PI-1500	\$43.25
5Fr x 8cm (3.25in)	with a 18Ga & 20Ga lumen, 5Fr x 1.5 inch peel away introducer	PI-1640	\$43.25
5Fr x 13cm (5.25in)	with a 18Ga & 20Ga lumen, 5Fr x 1.5 inch peel away introducer	PI-1600	\$43.25
5Fr x 25cm (10in)	with a 18Ga & 20Ga lumen, 5Fr x 2.5 inch peel away introducer	PI-1650	\$43.25

### Triple Lumen

5.5Fr x 8cm (3.25in)	with one 20Ga & two 22Ga lumens, 6Fr x 1.5 inch peel away introducer	PI-1740	\$46.30
5.5Fr x 13cm (5.25in)	with one 20Ga & two 22Ga lumens, 6Fr x 1.5 inch peel away introducer	PI-1700	\$46.30

PEEL AWAY INTRODUCER KIT  
includes peel away catheter-  
over-needle introducer, catheter,  
extension tube and tethered cap.



### Guidewire in Dispenser

One soft flexible 'J' end, one soft flexible straight end.

DESCRIPTION	ITEM #	PRICE
Sizes include:		
.018 x 45cm - Use with 4Fr & 5Fr catheters (MILA #1500, 1540, 1600, 1640, 1700, 1740, 1650, 1906, 1910, SA1606, SA1610)	GW1845	\$6.20
.032 x 45cm - Use with 16ga catheters (MILA #1606, 1606-2P, 1610-2P)	GW3245	\$6.20
.032 x 60cm - Use with 16ga catheter (MILA #1610, 1625-2P)	GW3260	\$6.20
.035 x 45cm - Use with 14ga and 7Fr catheters (MILA #1406, 1406-2P, 1410-2P, 1720)	GW3545	\$6.20
.035 x 60cm - Use with 14ga and 7Fr catheters (MILA #1410, 1620, 1630)	GW3560	\$6.20



### Peel Away Introducers


Facilitate the placement of catheters and tubing with these high quality introducers. They peel apart cleanly after introduction. Use in MILA long term catheter kits and PICC kits.

DESCRIPTION	ITEM #	PRICE
3Fr (20ga) x 3.75cm (1.5in)	4203	\$15.45
4Fr (18ga) x 3.75cm (1.5in)	4204	\$15.45
5Fr (16ga) x 3.75cm (1.5in)	4205.5	\$15.45
5Fr (16ga) x 6.25cm (2.5in)	4205	\$15.45
6Fr (14ga) x 3.75cm (1.5in)	4206.5	\$15.45
6Fr (14ga) x 6.25cm (2.5in)	4206	\$15.45



### Vessel Dilator



DESCRIPTION	ITEM #	PRICE
4Fr x 5cm	DIL4	\$2.50
5Fr x 6cm	DIL5	\$2.50
6Fr x 6cm	DIL6	\$2.50
7Fr x 10cm	DIL7	\$2.50
8Fr x 10cm	DIL8	\$2.50
9Fr x 10cm	DIL9	\$2.50



## Extended Use MILACATH®

- Kink resistant polyurethane catheter-over-needle.

DESCRIPTION	ITEM #	PRICE
20ga x 6.25cm (2.5in)	2011	\$10.20
18ga x 6.25cm (2.5in)	1811	\$10.20
16ga x 7.5cm (3.0in)	1603	\$10.20
16ga x 13cm (5.25in)	1611	\$10.95
14Ga x 9cm (3.5in)	1403	\$10.20
14Ga x 13cm (5.25in)	1411	\$10.95
12Ga x 13cm (5.25in)	1211	\$10.95



## Short Term Catheters

- Teflon catheter-over-needle.

DESCRIPTION	ITEM #	PRICE
22Ga x 45mm (1.75in)	2217	\$3.30
14Ga x 7cm (2.75in)	1402	\$4.60
12Ga x 13cm (5.25in)	1212	\$9.65
Fenestrated 12Ga x 13cm (5.25in) 4 side holes near distal end	1212F	\$10.80
10Ga x 7.5cm (3in)	1003	\$10.00
10Ga x 15cm (6in)	1006	\$10.50



## Introcan® Safety™ Catheters

- Introcan® Safety™ polyurethane catheters feature a safety clip that automatically covers the needle point when it is withdrawn, reducing accidental needle sticks.
- Push-off plate on hub for smooth insertion.
- Guarded flashback chamber.

DESCRIPTION	ITEM #	PRICE
24Ga x 19mm (.75in)	IC2475	\$1.95
22Ga x 25mm (1.0in)	IC2201	\$1.95
20Ga x 32mm (1.25in)	IC2001	\$1.95
18Ga x 45mm (1.75in)	IC1801	\$1.95



## Winged Introcan® Safety™ Catheters


DESCRIPTION	ITEM #	PRICE
24Ga x 19mm (.75in)	WIC2475	\$1.95
22Ga x 25mm (1.0in)	WIC2201	\$1.95
20Ga x 32mm (1.25in)	WIC2001	\$1.95
18Ga x 45mm (1.75in)	WIC1801	\$1.95



## Long Line Catheters

- Long term catheter may be inserted in any vein.
- Polyurethane kink-resistant catheter can be cut to desired length.
- Secure female luer lock is placed onto catheter after insertion into vessel.
- Moveable wing provides method of secure attachment.
- Kit includes catheter, introducers, moveable suture wing and female luer lock adapter.

DESCRIPTION	ITEM #	PRICE
20Ga catheter adjustable up to 45cm (18in)	LL2045	\$31.95
19Ga catheter adjustable up to 45cm (18in)	LL1945	\$18.25
19Ga catheter adjustable up to 90cm (36in)	LL1990	\$18.25
Replacement Clamshell for 20Ga catheter	8055	\$3.10
Replacement Clamshell for 19Ga catheter	8050	\$3.10



### Other Long Line Catheters

16Ga (5Fr) x 120cm (48in) catheter with peel away introducer.

See catalog page 7.

ITEM #/ PRICE	PI-1648 / \$20.85
---------------	-------------------

## Veincath Catheters

- Innovative replacement for the Venocath.™
- Catheter through needle insertion.
- **Protective pouch** to advance catheter with stylet.
- After catheter is placed, the needle is safely discarded before easy attachment of the female luer.
- Economical.

DESCRIPTION	ITEM #	PRICE
20Ga x 30cm (12in) catheter with 18ga introducer	V2030	\$13.50
Replacement Clamshell	8050	\$3.10



## PICC - Silicone Catheter Kit

- PICC - **P**eripherally **I**nserted **C**entral **C**atheter.
- May be inserted in any vessel.
- Silicone material for long term use of 30 + days.
- Catheter can be cut to desired length.
- Insertion method via peel away introducer.
- Catheter has a stylet wire to facilitate advancement.
- Kit includes catheter, peel away introducer and ruler.


DESCRIPTION	ITEM #	PRICE
2Fr (23ga) x 40cm (16in) Single Lumen	3240	\$56.45
3Fr (20ga) x 60cm (24in) Single Lumen	3360	\$56.45
4Fr (18ga) x 60cm (24in) Single Lumen	3460	\$56.45
5Fr (16ga) x 60cm (24in) Single Lumen	3560	\$56.45
5Fr x 60cm (24in) Double Lumen (2-18Ga)	3565	\$69.30



## CVP Manometer

35cm water manometer with four-way stopcock, 125cm (50 inch) extension, injection site and luer lock adapter.

ITEM # / PRICE	CVP35 / \$20.55
----------------	-----------------




## Extensions with Normally Closed Valve

### Features

- Catheter saving 'normally closed valve'.
- Bi-directional valve for infusion and sampling.
- High flow or small bore tubing.
- Spin luer lock.

DESCRIPTION	ITEM #	PRICE
High Flow - 13Ga x 45cm (18in)	8518	\$3.85
Small Bore T-Set - 16Gax17cm (7in)	8570	\$4.60
High Flow T-Set - 12Ga x 17cm (7in)	8572	\$4.60
Double High Flow T-Set - 12Ga	8575	\$4.60
Small Bore - 16Ga x 17cm (7in)	8585	\$3.85
High Flow - 13Ga x 17cm (7in)	8590	\$3.85
Slide Clamp for 8585 & 8590 (50 per bag)	8599	\$5.90




### Specialty Extension Sets

DESCRIPTION	ITEM #	PRICE
Pressure Monitoring Extension 16Ga x 180cm (72in) Includes female luer lock and male spin luer lock.	8515	\$3.75
Large Bore High Flow 10Ga x 180cm (72in) Includes female luer lock, injection port and male spin luer lock.	8520	\$6.50
Small Bore Extension 20Ga x 100cm (40in) Includes female and male luer lock.	8540	\$3.85

## Extensions with MicroCLAVE®

- One-piece, closed, needle-free connector.
- Approved for saline flush.
- Heparin not required.

DESCRIPTION	ITEM #	PRICE
Smallbore T-Set with MicroClave, Injection Site, Clamp, Rotating Luer Lock (18cm/7in length)	MC40	\$4.85
Smallbore T-Set with dual MicroClaves, Injection Site, Clamp, Rotating Luer Lock (18cm/7in length)	MC43	\$5.00
Smallbore Bifuse 'Y' Set with three Clamps, Rotating Luer Lock (18cm/7in length)	MC44	\$5.00
MicroClave with Luer Lock	8098	\$2.00
Stemless Luer Lock Caps (10 pack)	8099	\$3.75



## Caps

- **Safety Injection Cap** (8080) has a built-in stop to prevent accidental tubing or technician puncture.
- **Tethered Cap** (8090) attaches to the catheter or extension hub to cover the injection site.
- **Needleless Injection Cap** (8095 & 8090) for direct attachment of an IV administration line or syringe. Includes a tethered cap (8090).
- **MicroClave** (8098) is a one-piece, closed, needle-free connector. Labeled for saline flush. Heparin not required.
- **Stemless Luer Lock Cap** (8099) 10 pack, non-sterile.
- **Clamshell Female Luers** (8050, 8055)  
Hinged female luer lock for attaching after catheter placement.

DESCRIPTION	ITEM #	PRICE
Safety Injection Cap	8080	\$0.80
Tethered Cap	8090	\$0.80
Needleless Injection Cap	8095	\$3.00
MicroClave	8098	\$2.00
Stemless Luer Lock Cap	8099	\$3.75
Clamshell - 18Ga Female Luer	8050	\$3.10
Clamshell - 20Ga Female Luer	8055	\$3.10




## Instructional DVD

Demonstrates placement of:

- Centesis Valve
- Diffusion/Wound Catheter
- Endotracheal Tubes - lateral and dorsal
- Epidural Catheter
- Foley Catheter
- Infu-Disk
- Nasogastric Feeding Tube
- Oxygen Catheter - nasopharyngeal & nasotracheal
- Peritoneal Lavage
- Thoracentesis
- Tracheostomy Tube
- Venous Access Catheters:
  - Guidewire insertion technique
  - Peel away insertion technique
  - Over-the-wire Catheter Exchange
  - PIC catheter

ITEM # / PRICE

DVD100 / \$20.00

**MILA**   
INTERNATIONAL, INC.  
*Medical Instrumentation for Animals*  
859.371.1722 \* Toll Free 888.MILA.INT \* www.milaint.com

**INSTRUCTIONAL DVD**  
**REORDER NUMBER: DVD100**

**Demonstrates placement of:**


- Centesis Valve
- Diffusion/Wound Catheter with Infu-Disk
- Endotracheal Tube - lateral and dorsal
- Epidural Catheter
- Feeding Tube - nasoesophageal & nasogastric
- Foley Catheter
- Infu-Disk
- Oxygen Catheter - nasopharyngeal & nasotracheal
- Peritoneal Lavage
- Thoracentesis
- Tracheostomy Tube using Endotracheal Tube
- Venous Access Catheters:
  - Long term - guidewire insertion and peel away insertion
  - Catheter Exchange - Over-the-wire technique
  - Peripherally Inserted Central Catheter (PICC)

## Tomcat / Small Animal Urinary Catheter


- Same catheter used to unblock and leave indwelling.
- Improvement to the standard Tomcat, Jackson and Slippery Sam™ urinary catheters.
- Non-irritating polyurethane material.
- Rigid for easy insertion.
- Softens indwelling.
- Finished tip.
- Luer lock.
- Moveable suture wing allows length adjustment.
- Optional stylet.

DESCRIPTION	ITEM #	PRICE
3.5Fr - adjustable up to 15cm (6in)	UC306	\$10.25
3.5Fr - adjustable, with stylet	UC306S	\$12.90
3.5Fr - adjustable up to 25cm (10in)	UC310	\$10.25
3.5Fr - adjustable, with stylet	UC310S	\$12.90
5Fr - adjustable up to 15cm (6in)	UC506	\$12.90
5Fr - adjustable, with stylet	UC506S	\$15.30
5Fr - adjustable up to 25cm (10in)	UC510	\$12.90
5Fr - adjustable, with stylet	UC510S	\$15.30

See page 15 for Collection Bag



Tomcat catheters with moveable wing and optional stylet



Can be ordered with or without a stylet



## Foley Urinary Catheters

Radiopaque silicone foley catheter with balloon cuff and wire stylet. Water activated lubricant on 6 French and 8 French sizes.

DESCRIPTION	ITEM #	PRICE
6Fr x 60cm (24in) with 1.5cc balloon, suture wing	UC624	\$37.00
6Fr x 90cm (36in) with 1.5cc balloon, suture wing	UC636	\$37.00
8Fr x 30cm (12in) with 3cc balloon, suture wing	UC812	\$37.00
8Fr x 60cm (24in) with 3cc balloon, suture wing	UC824	\$37.00
8Fr x 90cm (36in) with 3cc balloon, suture wing	UC836	\$37.00
10Fr x 45cm (18in) with 5cc balloon	UC1018	\$37.00
10Fr x 90cm (36in) with 5cc balloon	UC1036	\$37.00
12Fr x 30cm (12in) with 10cc balloon	UC1212	\$40.00
12Fr x 55cm (21in) with 10cc balloon	UC1221	\$40.00
14Fr x 23cm (9in) with 10cc balloon	UC149	\$40.00
16Fr x 55cm (21in) with 20cc balloon	UC1621	\$46.00
24Fr x 42cm (16in) with 30cc balloon	UC2416	\$52.00
24Fr x 150cm (59in) with 30cc balloon	UC2459	\$52.00



## Urine Collection Systems

DESCRIPTION	ITEM #	PRICE
1000cc Collection Bag with male luer lock and twist drain.	UC1000	\$18.00
Urine Collection Kit - Closed sterile system collection bag with sampling port, 72in collection line, injection port and female Christmas tree adapter.	UC1000K	\$23.00



## Ferret Urinary Catheter

- Two Styles - Closed end with side ports.  
Open ended with side ports.
- Non-irritating polyurethane material.
- Length adjustable by cutting the proximal end.
- Attach female luer lock after catheter placed.
- Suture wing.


DESCRIPTION	ITEM #	PRICE
3Fr x 20cm (8in) closed end	UCC308	\$31.00
3Fr x 20cm (8in) open end	UCO308	\$31.00



**Nasogastric Feeding Tubes**

- Marked, non-irritating radiopaque polyurethane.
- All tubes have a stylet.
- Y-adapter for feeding, medication and flushing.
- Suitable for oropharyngeal oxygen.

DESCRIPTION	ITEM #	PRICE
5Fr x 55cm (22in)	NG522S	\$19.45
6Fr x 55cm (22in)	NG622	\$15.45
6Fr x 90cm (36in)	NG636	\$15.45
8Fr x 108cm (43in)	NG843	\$15.45
8Fr x 140cm (55in)	NG855	\$15.45
10Fr x 90cm (36in)	NG1036	\$15.45
10Fr x 108cm (43in)	NG1043	\$18.00
10Fr x 140cm (55in)	NG1055	\$18.00
12Fr x 108cm (43in)	NG1243	\$19.55
14Fr x 125cm (50in)	NG1450	\$29.70
Length adjustable and includes nose piece		
18Fr x 250cm (100in)	NG18100	\$32.35
Length adjustable and includes nose piece		
24Fr x 300cm (120in)	NG24120	\$36.00
For enteral, extruded and other feeds soaked and made into a fine slurry		
24Fr x 240cm (96in)	NG2496	\$35.75
Hard teflon material, no stylet		



**Weighted Nasogastric Feeding Tubes**

- Marked, non-irritating radiopaque polyurethane.
- All tubes have a stylet.
- Y-adapter for feeding, medication and flushing.

DESCRIPTION	ITEM #	PRICE
6Fr x 55cm (22in) - 1 Gram	NG622W	19.55
8Fr x 108cm (43in) - 3 Gram	NG843W-3	21.85
8Fr x 108cm (43in) - 7 Gram	NG843W-7	21.85
10Fr x 90cm (36in) - 3 Gram	NG1036W	26.25
12Fr x 108cm (43in) - 3 Gram	NG1243W	26.55



## Esophagostomy Feeding Tubes

### POLYURETHANE MATERIAL

- Non-irritating radiopaque polyurethane material.
- Y-adapter for feeding, medication and flushing.
- Length adjustable **after** placement by cutting the proximal end.
- Secure Y-adapter.

DESCRIPTION	ITEM #	PRICE
Feline - 14Fr adjustable to 33cm(13in)	E1413	\$20.40
Canine - 14Fr adjustable to 75cm (30in)	E1430	\$20.40
Small Animal - 14Fr adjustable to 125cm (50in), includes stylet	E1450	\$22.65
Small Animal - 18Fr adjustable to 75cm (30in)	E1830	\$22.65
Y-Repair Kit	E-Y	\$9.00



## Esophagostomy Feeding Tubes

### SILICONE MATERIAL

- Radiopaque material.
- Open end with multiple side holes.
- Female luer lock adapter.

DESCRIPTION	ITEM #	PRICE
12Fr x 38cm (15in) with 5cm of sideholes	SE1215	\$21.00
14Fr x 38cm (15in) with 7cm of sideholes	SE1415	\$21.00
20Fr x 38cm (15in) with 10cm of sideholes	SE2015	\$21.00
20Fr x 55cm (22in) with 10cm of sideholes	SE2022	\$21.00
30Fr x 38cm (15in) with 12cm of sideholes	SE3015	\$21.00
30Fr X 55cm (22in) with 12cm of sideholes	SE3022	\$21.00



## Cubby – Low Profile Gastrostomy Device

- Balloon style cuff which enables insertion and removal without trauma after a stoma is formed.
- Can be easily replaced by the client.
- Kit includes patient care kit with a 6cc & 60cc syringes, right angle feeding set, bolus feeding set and care giver guide.
- Available in 14Fr, 16Fr, 18Fr, 20Fr and 24Fr sizes and lengths of 1.5cm, 2.0cm, 2.5cm, 3.0cm, 3.5cm, 4.0cm and 4.5cm.
- \$133.90 each.
- Stoma measuring device available (#M5555)/ \$9.80 each.



## Gastrostomy Kits – Endoscope or Surgically Placed


**Surgical Kit** includes 30 inch peg tube, right angle fixation device and Y-adapter.

DESCRIPTION	ITEM #	PRICE
16Fr Surgically Placed	PGS16	\$48.50
20Fr Surgically Placed	PGS20	\$48.50

### PEG - Percutaneous Endoscopic Gastrostomy

- Radiopaque polyurethane material is non-irritating and stronger than silicone.
- Collapsible bumper to facilitate insertion and removal. Less traumatic removal than mushroom style.
- Right angle fixation device that holds tube securely.
- Locking Y-adapter for feeding, irrigation, decompression and medication without disconnection.
- Not recommended for large dogs.
- Kit includes introduction needle, guidewire, 36 inch peg tube, right angle fixation device and Y-adapter.

DESCRIPTION	ITEM #	PRICE
16Fr Endoscopically Placed	PEG16	\$68.50
20Fr Endoscopically Placed	PEG20	\$68.50
Fixation Bar - 16Fr	FIX16	\$3.30
Fixation Bar - 20Fr	FIX20	\$3.30
Right Angle Fixation Device - 16Fr	RAFIX16	\$6.85
Right Angle Fixation Device - 20Fr	RAFIX20	\$6.85
Y-Repair Kit - 16Fr	RPEG16	\$16.05
Y-Repair Kit - 20Fr	RPEG20	\$16.05




1. PEG Y  
2. PEG Bumper  
3. PEG Tube  
4. PEG Wire  
5. PEG Assembly

## Jejunostomy Tubes

- Must be placed through either a 16Fr or 20Fr PEG tube.
- Continuous nutritional support via a jejunostomy tube positioned in duodenum or jejunum.
- Silicone plug and capped stylet prevent the escape of air during placement.
- Sealed Twoomey adapter permits evacuation of stomach from side port and placement of J-tube through top port.

DESCRIPTION	ITEM #	PRICE
6Fr x 36in (Order PEG16 or PGS16)	J636	\$56.65
8Fr x 36in (Order PEG20 or PGS20)	J836	\$56.65
10Fr x 36in (Order PEG20 or PGS20)	J1036	\$56.65




## Needles – Spinal & Tuohy

Spinal and tuohy needles with marked and notched hub for easy tip orientation.

DESCRIPTION	ITEM #	PRICE
<b>Spinal Needles</b>		
22Ga x 15cm (6in)	2206	\$8.05
20Ga x 20cm (8in)	2008	\$9.25
18Ga x 15cm (6in)	1806	\$7.10
18Ga x 20cm (8in)	1808	\$9.80
18Ga x 25cm (10in)	1810	\$10.30
<b>Tuohy Needles</b>		
18Ga x 7.5cm (3in)	T1803	\$9.25
18Ga x 15cm (6in)	T1806	\$10.50
20Ga x 6.25cm (2.5in)	T2025	\$9.25
20Ga x 8.75cm (3.5in)	T2035	\$9.25



Spinal Needles



Tuohy Needles

## Bone Marrow Aspiration Needle

11Ga x 10cm (4in) double diamond aspiration needle, obturator pushing rod, probe guide and luer lock cap.

ITEM # / PRICE	BMAN / \$27.80
----------------	----------------



## Controlled Flushing Set

- Set includes vented bag spike, 90cm (36in) tubing and bi-directional valve.
- Use for various flushing needs - joint, ear, wound, ovum pickup.
- The bi-directional valve allows for repeat controlled flushing and vacuuming without opening the system.
- Eliminates the need to pre-fill syringes.
- Spike the flush medium, attach syringe to female luer of valve, fill syringe, flush, then repeat process.

ITEM # / PRICE	CFS36 / \$14.50
----------------	-----------------




## Guidewire Inserted Chest Tube

- Single lumen polyurethane catheter with multiple fenestrations.
- Use for pyothorax, hemothorax, pneumothorax, emergency peritoneal lavage and pericardial drain.
- Insertion using the catheter-over-guidewire technique.
- Normally closed valve with luer adapters prevents air from entering chest.
- Additional suture wing for secure adjustment and attachment.
- Multiple drainage holes. 14Ga fenestrated up to 8cm mark, 12Ga fenestrated up to 17cm mark.
- Kits include introducers, 60cm guidewire, catheter, suture wings, normally closed valve and tethered cap.

DESCRIPTION	ITEM #	PRICE
14Ga x 20cm (8in)	CT1410	\$29.00
12Ga x 30cm (12in)	CT1210	\$32.00



Photos courtesy of Mike Ethier DVM, DVSc, DACVECC




See MILA website for discussion by Mike Ethier DVM

## Silicone Chest Tubes

- Single lumen open-end with multiple side holes.
- Use for drainage of air or fluid.
- Surgically placed.
- Includes clamp, barbed adapter or stopcock.

DESCRIPTION	ITEM #	PRICE
12Fr x 38cm (15in) with 5cm of sideholes	CT1215	\$18.45
14Fr x 38cm (15in) with 7cm of sideholes	CT1415	\$18.45
20Fr x 38cm (15in) with 10cm of sideholes	CT2015	\$21.85
20Fr x 55cm (22in) with 10cm of sideholes	CT2022	\$21.85
30Fr x 38cm (15in) with 12cm of sideholes	CT3015	\$21.85
30Fr x 55cm (22in) with 12cm of sideholes	CT3022	\$21.85

See page 34 for barbed adapters



## Chest Tube Suction Bulb

- Continuous gentle suction.
- Luer lock adapter for secure attachment to chest tube.
- Volume marking provides easy measurement.
- Non-reflux valve.

DESCRIPTION	ITEM #	PRICE
100cc suction bulb	CG100	\$18.00
150cc suction bulb	CG150	\$18.00
200cc suction bulb	CG200	\$18.00
400cc suction bulb	CG400	\$22.00






## Chest Drain Adapters

### CHEST DRAIN/CENTESIS ADAPTER

**Eliminates the need for stopcocks!** Unique valve allows aspiration and evacuation with a syringe. Luer lock adapter for a needle. Includes a barbed adapter for 18-34Fr chest tubes. Also available with an in-line injection site.


DESCRIPTION	ITEM #	PRICE
Spin Luer lock only - no tubing	175440	\$9.80
Total Length 12.5cm (5in)	175450	\$13.00
Total Length 40cm (16in)	175460	\$14.15
In-Line Injection Site (10in)	175470	\$14.85



### CHEST TUBE ADAPTER with stopcock

Adapter assembled to include a three-way stopcock, clamp and barbed adapter for 18-34Fr tubes.

DESCRIPTION	ITEM #	PRICE
16Fr x 17cm (7in)	175420	\$9.55
20Fr x 25cm (10in)	175425	\$9.55



## 'Heimlich' Valve

One way valve with barbed adapters. Small, five inches long and lightweight.


ITEM # / PRICE	H2900 / \$7.65
----------------	----------------



## Suction Connector for Bilateral Chest Tubes

This set connects bilateral chest tubes to a single vacuum source. Luers allow easy disconnect and reconnect.

DESCRIPTION	ITEM #	PRICE
Small Diameter (20Fr) Suction Set	175800	\$13.55
Small Diameter Suction Kit	175810	\$28.50
Includes two 175420 adapters		
Large Diameter (30Fr) Suction Set	175820	\$13.55
Large Diameter Suction Kit	175830	\$28.50
Includes two 175425 adapters		



**For demonstration, order DVD100**

## Grenade – Continuous Suction

### Multipurpose

- Two options - Introducer, drain and adapter or Introducer, drain, adapter and 100cc suction bulb.
- 90cm (36in) round silicone material.
- Radiopaque line.
- Multiple perforations over 25cm (10in) of distal end.


DESCRIPTION	ITEM #	PRICE
8Fr drain with introducer	D890	\$12.35
10Fr drain with introducer	D1090	\$12.35
8Fr drain, introducer, adapter and 100cc suction bulb	DK890	\$27.35
10Fr drain, introducer, adapter and 100cc suction bulb	DK1090	\$27.35



### Jackson Pratt

Surgically placed silicone drain with full perforations.


DESCRIPTION	ITEM #	PRICE
4.5mm	JP4	\$10.00
7mm	JP7	\$10.00
10mm	JP10	\$10.00
13mm	JP13	\$10.00



### Blake

Full fluted. Ideal for peritoneal dialysis.

DESCRIPTION	ITEM #	PRICE
8mm Flat	BL8	\$15.00
11mm Flat	BL11	\$15.00



- Low pressure design provides continuous gentle suction and minimizes tissue trauma.
- Volume marking provides easy measurement.
- Non-reflux valve.
- Flat bottom design for easy positioning and clip for secure attachment.

DESCRIPTION	ITEM #	PRICE
100cc suction bulb w/one-way valve	G100	\$15.00
150cc suction bulb w/one-way valve	G150	\$15.00
200cc suction bulb w/one-way valve	G200	\$15.00
400cc suction bulb w/one-way valve	G400	\$20.00




See page 20 for Chest Tube Suction Bulb with leur adapter.

## BAL, TTL and Tracheal Wash Kit

### Small Animal Trans-Tracheal Wash Kit

Catheter-over-needle insertion method. The semi-rigid flushing catheter has a stylet and luer adapter.

DESCRIPTION	ITEM #	PRICE
<b>Cats</b> - 19Ga x 15cm (6in) flushing catheter with 16Ga x 1.25in teflon introduction catheter	TW1900	18.00
<b>Small Dog and Cats</b> - 16Ga x 25cm (10in) flushing catheter with 14Ga x 2.75in teflon introduction catheter	TW1610	18.00
<b>Medium Sized Dog</b> - 16Ga x 41cm (16in) flushing catheter with 14Ga x 2.75in teflon introduction catheter	TW1616	23.55
<b>Catheter Only</b> - 16Ga x 41cm (16in) flushing catheter <i>without stylet</i> with 14Ga x 2.75in teflon introduction catheter	TW1617	17.00
<b>Large Dog</b> - 16Ga x 70cm (28in) flushing catheter with 14Ga x 2.75in teflon introduction catheter	TW1628	23.55



INTRODUCERS are catheter-over-needle


*For Large Animal, see page 32*

## Endoscope Catheters

Use for BAL, TTL and Tracheal Wash


SINGLE LUMEN catheter with female luer lock.

DESCRIPTION	ITEM #	PRICE
2.0mm x 60cm (24in)	EDC1624	\$10.75
2.0mm x 120cm (48in)	EDC1648	\$10.75



DOUBLE LUMEN enables sterile specimen retrieval. Outer catheter has a glycol plug which is dispelled when advancing the inner catheter for flush and specimen retrieval.

DESCRIPTION	ITEM #	PRICE
2.0mm x 90cm (36in)	EMAC090	28.05




*For equine BAL, see page 32*

## Canine GDV / Bloat Tube

- Dual lumen tube for simultaneous flushing and evacuation of the stomach especially during surgery for food bloat.
- Four side holes.
- Heavy duty silicone material.


DESCRIPTION	ITEM #	PRICE
Giant Breed 30mm (1.2in) OD x 180cm (72in)	GDV1	\$139.05
Large Breed 25mm (1.0in) OD x 180cm (72in)	GDV2	\$139.05
Medium Breed 17mm (0.6in) OD x 125cm (50in)	GDV3	\$110.00



## Epidural Pain Management Kit

- Cut to desired length – catheter will not crush.
- Long term pain management.
- Radiopaque multi-port catheter (90cm/36in) with stylet.
- Use with MILA administration pumps (see page 26/27).
- Kit includes catheter, tuohy introducer needle, loss of resistance syringe, filter, luer lock catheter clamp, suture wing, blade and drape.

DESCRIPTION	ITEM #	PRICE
20Ga catheter - 18Ga x 7.5cm (3in) tuohy needle	9203	\$38.45
20Ga catheter - 18Ga x 15cm (6in) tuohy needle	9206	\$38.45




See MILA website for detailed insertion instructions written by Shane Bateman DVM, DVSc, DACVECC.

## Diffusion Catheter / Wound Catheter

- 16Ga x 30cm (12in) polyurethane catheter with micro ports in the distal end to evenly disperse fluid over 2 inches, 4 inches, 6 inches, 7.5 inches or 9 inches.
- For delivery of pain medication or antibiotics directly to a wound.
- Can be surgically placed at the wound site or by using a peel away introducer.
- Use with MILA administration pumps (see page 26/27), a syringe, mini pump, or IV pump for continuous infusion.

DESCRIPTION	ITEM #	PRICE
2in (5cm) Dispersion Area Catheter length = 11cm (4.5in)	DIFF2	\$20.00
4in (10cm) Dispersion Area	DIFF4	\$20.00
6in (15cm) Dispersion Area	DIFF6	\$20.00
7.5in (18cm) Dispersion Area	DIFF7	\$20.00
9in (22.5cm) Dispersion Area	DIFF9	\$20.00
.2 Micron Filter (10 pack) Individually packed / sterile	2420	\$47.00



See MILA website for sample protocols by Emily McCobb DVM MS DACVA



Photo courtesy of William Muir DVM, PhD, The Ohio State University



## Non-Cuffed

- MRI Compatible
- Reuseable
- Ethylene Oxide Sterilization Recommended

Proximal fitting accepts a 15mm adapter.

DESCRIPTION	ITEM #	PRICE
1.0mm ID x 9.5cm tube with wire stylet	ETN10	\$14.05
1.5mm ID x 10cm tube with wire stylet	ETN15	\$14.05



## Stepped

- MRI Compatible
- Reuseable
- Ethylene Oxide Sterilization Recommended

DESCRIPTION	ITEM #	PRICE
2.0mm ID, 5.0mm OD - 9Fr x 13cm	ETN20	\$14.05
2.5mm ID, 6.3mm OD - 12 Fr x 13cm	ETN25	\$14.05
3.0mm ID, 8.0mm OD - 13.6Fr x 15cm	ETN30	\$14.05
3.5mm ID, 8.6mm OD - 16Fr x 15cm	ETN35	\$14.05
4.0mm ID, 12mm OD - 18.3Fr x 15cm	ETN40	\$14.05



## Cuffed

- MRI Compatible
- Reuseable
- Ethylene Oxide Sterilization Recommended

DESCRIPTION	ITEM #	PRICE
2mm ID, 4.0mm OD - 12Fr x 14cm (5.5in)	ETC20	\$41.20
3mm ID, 4.7mm OD - 14Fr x 16cm (6.4in)	ETC30	\$41.20
4mm ID, 6.0mm OD - 18Fr x 18cm (7.2in)	ETC40	\$41.20
5mm ID, 7.3mm OD - 20Fr x 20cm (8in)	ETC50	\$41.20
6mm ID, 8.7mm OD - 28Fr x 28cm (11.2in)	ETC60	\$41.20
7mm ID, 10mm OD - 30Fr x 55cm (22in)	ETC755	\$44.30
8mm ID, 11mm OD - 33Fr x 55cm (22in)	ETC855	\$44.30
9mm ID, 12mm OD - 33Fr x 33cm (13in)	ETC933	\$44.30
9mm ID, 12mm OD - 37Fr x 55cm (22in)	ETC955	\$62.85
10mm ID, 14mm OD - 43Fr x 33cm (13in)	ETC1033	\$44.30
10mm ID, 14mm OD - 43Fr x 55cm (22in)	ETC1055	\$62.85
11mm ID, 15mm OD - 47Fr x 33cm (13in)	ETC1133	\$44.30
11mm ID, 15mm OD - 47Fr x 55cm (22in)	ETC1155	\$64.90
12mm ID, 17mm OD - 51Fr x 33cm (13in)	ETC1233	\$44.30
12mm ID, 17mm OD - 51Fr x 55cm (22in)	ETC1255	\$64.90



## Spring Powered Syringe Pump

- Economical.
- Spring driven - No battery or electrical requirements.
- Saves time to administer proper dosage or multiple dosing.
- Decreased use of high-tech capital equipment for non high-tech situations.
- Lower cost per use.
- Easy refill without disconnect.
- Accuracy rate +/- 15 %.

DESCRIPTION	ITEM #	PRICE
Syringe Pump – 10ml	SF10	\$41.20
Refill Port with Extension Set	FCT REV	\$9.55
10ml winged syringe	SYR10	\$1.35

### Flow Control Tubing

10ml in 15 minutes – 40ml/hour	FCT15	\$9.55
10ml in 30 minutes – 20ml/hour	FCT30	\$9.55
10ml in 2 hours – 5ml/hour	FCT2	\$9.55
10ml in 4 hours – 2.5ml/hour	FCT4	\$9.55
10ml in 12 hours – 0.8ml/hour	FCT12	\$11.55



Courtesy of D. Mader, MS, DVM, DAVPB, Marathon Veterinary Hospital

## Elastomeric Pump

- Four layer balloon pump.
- Recoil of balloon forces fluid from pump.
- 100ml fluid reservoir.
- Luer lock flow control tube, 0.8ml/hr, is attached to pump.
- Optional FCT4 (10ml/4hrs) included in package. 0.8ml/hr tube must be removed from pump before attaching FCT4.
- Flow rate can be changed by using other flow control tubes.
- Accuracy rate +/- 15 %.

ITEM # / PRICE	7100 / \$37.10
----------------	----------------

For other flow rates, replace 0.8ml/hr tube with one of the following

DESCRIPTION	ITEM #	PRICE
10ml in 15 minutes – 40ml/hour	FCT15	\$9.55
10ml in 30 minutes – 20ml/hour	FCT30	\$9.55
10ml in 2 hours – 5ml/hour	FCT2	\$9.55
10ml in 4 hours – 2.5ml/hour	FCT4	\$9.55
10ml in 12 hours – 0.8ml/hour	FCT12	\$11.55
.2 micron filter (10 pack)	2420	\$46.35



## Infu-Disk™

- 10ml fluid reservoir.
- Choose desired delivery time.
- Single use and disposable.
- Lightweight - 19 grams.
- Flow rate not affected by fluid viscosity.
- Large bore outlet tubing prevents occlusion.
- Battery powered - fill the fluid reservoir and press the start button.
- Accuracy rate better than +/- 10%.




DESCRIPTION	ITEM #	PRICE
0.06ml/hr – 7 day delivery	ID06	\$57.00
0.08ml/hr – 5 day delivery	ID08	\$57.00
0.10ml/hr – 4 day delivery	ID010	\$57.00
0.14ml/hr – 3 day delivery	ID14	\$57.00
0.21ml/hr – 2 day delivery	ID21	\$57.00
0.42ml/hr – 1 day delivery	ID42	\$57.00
0.50ml/hr – 20 hour delivery	ID50	\$57.00
1.0ml/hr – 10 hour delivery	ID10	\$57.00
T-Adapter	ID-T	\$3.20
W-Adapter	ID-W	\$3.20

**For demonstration, order DVD100**

## Nasolacrimal Duct Catheter

Indwelling catheter provides flushing and medication delivery. Catheters include a stylet to improve rigidity during placement.

DESCRIPTION	ITEM #	PRICE
Open end		
3.5Fr (19Ga) x 25cm (10in)	NL325	\$15.50
5Fr (16Ga) x 25cm (10in)	NL525	\$15.50
5Fr (16Ga) x 25cm (10in) w/flexible stylet	NL525FS	\$15.50
Closed end, bullet tip with large side hole		
5Fr (16Ga) x 55cm (22in)	NL555	\$19.45




## Eye Lavage Kit / SPL

- For upper or lower lid.
- Single stick insertion method.
- Soft non-irritating disk and tubing.
- Allows application of eye medication away from the site.
- Easily used in clinical and farm setting.
- Kit includes 12Ga straight or 10° angled trocar, catheter, eyelid lifter, suture devices, extension tube, luer lock adapter and tethered cap.
- Use the elastomeric administration pump or Infu-Disk™ for continuous infusion.

DESCRIPTION	ITEM #	PRICE
Subpalpebral Eye Lavage Kit		
36in catheter	6612	\$58.00
Subpalpebral Eye Lavage Kit		
60in catheter	6612L	\$58.00
Subpalpebral Eye Lavage Kit		
Angled 10° trocar, 36in catheter	6613	\$58.00
Subpalpebral Eye Lavage Kit		
Angled 10° trocar, 60in catheter	6613L	\$58.00
Repair Kit	6615	\$13.00

Infu-Disk™ Pump recommended for continuous infusion




## Tracheal Wash / Aspiration Catheters

### Single Lumen Catheters

Single lumen catheter with female luer lock for medication delivery and sampling at any endoscope accessible site.


DESCRIPTION	ITEM #	PRICE
2.5mm x 190cm (76in)	EDC190	\$10.75
2.8mm X 400cm (158in)	EDC400	\$61.80
Heavier wall with two side ports		




### Double / Triple Lumen Catheters

Two and three lumen sampling catheters enable a sterile specimen retrieval via an endoscope. To maintain sterility during placement, these sets contain a glycol plug in the outer catheter. Once placed, the inner catheter is advanced for retrieval of the sterile specimen.


DESCRIPTION	ITEM #	PRICE
Two stage catheter - 2.0mm x 190cm (76in)	EMAC700	\$28.05
Two stage catheter - 2.0mm x 250cm (100in)	EMAC700L	\$33.65
Triple stage catheter - 2.5mm x 190cm (76in)	EMAC800	\$30.55



Double Lumen




Triple Lumen



## Lance-A-Lot

Aspiration and injection needle (17ga) with female luer lock. Ideal for lancing and injecting ethmoid hematomas and uterine cysts.

DESCRIPTION	ITEM #	PRICE
2.5mm x 190cm (76in)	LAN190	\$36.40



## Low Dose Insemination Set

Two stage catheter used for low dose insemination via an endoscope. The outer catheter is 2.0mm x 190cm and the 2 French inner catheter is 200cm in length. Includes 3cc luer lock syringe.

DESCRIPTION	ITEM #	PRICE
2.0mm x 200cm (80in)	INS200	\$46.35




## Indwelling Guttural Pouch Catheter Kit

Placed by viewing the guttural pouch through an endoscope up one nostril and placing the catheter introducer up the other nostril. Non-irritating polyurethane catheter for extended treatment of the guttural pouch. Kit includes polyurethane catheter with coils, bendable introducer and attachable luer lock adapter. The catheter can be left in for days, secured to the halter and protected by a grazing muzzle or cut to desired length.

Catheter dimension: 12Ga x 135cm (54 Inches)

DESCRIPTION	ITEM #	PRICE
Catheter Kit	9500	\$61.10
Catheter Only	9510	\$46.35
Introducer Only	9520	\$20.10




## Uterine Flush

### EMBRYO TRANSFER CATHETER

- Open end silicone catheter.
- Soft 50cc balloon cuff ideal for embryo transfer.
- Y-adapter tube assembly includes line to fluid source, drain line and graduated adapter line to catheter.
- Reuseable.


DESCRIPTION	ITEM #	PRICE
34Fr x 80cm (32in) catheter with 50cc balloon cuff	UFC3480	\$79.30
40Fr x 80cm (32in) catheter with 50cc balloon cuff	UFC4080	\$79.30
34Fr x 110cm (44in) cleansing flush catheter without balloon cuff	CF34110	\$50.45
Y-adapter tube assembly	UFYA	\$13.55
High Flow Fluid Transfer - 30cm (see pg 24)	175490	\$8.10
High Flow Fluid Transfer - 90cm (see pg 24)	175495	\$8.60



## Oocyte Collection Set

- For equine transvaginal oocyte collection.
- Laminar flow, no ridges on collection side.
- Aspiration needle is 12Ga x 60cm (24in), stainless steel, double lumen, with an ultra sharp etched tip to improve visualization under ultrasound.
- Flushing and collection lines are 14Ga x 180cm (72in) pvc attached to a 'Y' adapter.
- #6 silicone bottle stopper with collection port and vent port.

DESCRIPTION	ITEM #	PRICE
Oocyte Collection Set - Etched Tip Includes needle assembly, flushing and collection lines and bottle stopper assembly	OCS	\$75.20
Outer Replacement Needle - Etched Tip 12Gauge x 60cm	OCRN	\$29.85
Oocyte Collection Set - Non -Etched Tip Includes needle assembly, flushing and collection lines and bottle stopper assembly	OCS-2	\$75.20
Outer Replacement Needle - Non-Etched Tip 12Gauge x 60cm	OCRN-2	\$29.85
Oocyte Collection Needle Assembly 16Ga x 62.5cm (25in) needle with flushing and collection lines and bottle stopper assembly	OCN16	\$43.25
Replacement Lines Includes Y-port with attached flushing and collection lines	OCRL	\$8.75
Stopper Assembly Includes #6 silicone stopper with collection and vent ports	OCSA	\$22.65



## Low Dose Insemination Set

Two stage catheter used for low dose insemination via an endoscope. The outer catheter is 2.0mm x 190cm and the 2 French inner catheter is 200cm in length. Includes 3cc luer lock syringe.


DESCRIPTION	ITEM #	PRICE
2.0mm x 200cm (80in)	INS200	\$46.35




## Fluid Collection and Delivery

### High Flow IV Administration Sets


DESCRIPTION	ITEM #	PRICE
High flow set includes drip chamber, pinch clamp and spin luer lock. Total length is 72.5cm (29in).	7350	\$6.80




DESCRIPTION	ITEM #	PRICE
High flow set with vented spike, pinch clamp and spin luer lock. Use with fluid bags and/or bottles. Total length is 72.5cm (29in).	7360	\$6.80



DESCRIPTION	ITEM #	PRICE
High flow set includes drip chamber, six foot (180cm) tubing, pinch clamp, injection site and spin luer lock.	175497	\$7.85




DESCRIPTION	ITEM #	PRICE
High flow infusion set includes vented spike, six foot (180cm) tubing, injection site, clamp and spin luer lock.	175498	\$6.85



### EquiSet IV Administration Set

DESCRIPTION	ITEM #	PRICE
<i>Improved - Now available with injection site</i> Economical large bore (30Fr x 300cm/120in) non-vented Y-set with 2 spikes, 3 clamps, roller flow adjustment and spin luer lock adapter.	8V0002	\$10.20




### Secondary Venoset

DESCRIPTION	ITEM #	PRICE
Deliver a small amount of medication quickly. This set allows pressurization of a bag or bottle. Replacement for the discontinued Abbott secondary venoset. Tubing length is 23in (58cm).	175720	\$4.40




### Fluid Control Dial

DESCRIPTION	ITEM #	PRICE
Easily adjust infusion rate from 10ml/hour up to 250ml/hour when used with gravity system. When used with pressurized system, flow is four times faster. Tube length is 18in (45cm).	175480	\$8.90



### Blood Collection Set

DESCRIPTION	ITEM #	PRICE
High flow collection set that attaches directly to venous catheter and vacutainer bottle. Set includes six foot (180cm) tubing, injection site, chamber and spin luer lock.	175484	\$9.70



## High Flow Fluid Transfer Set


Ideal for uterine flush. Includes high flow spike, clamp and graduated adapter for 20-40Fr tubes.

DESCRIPTION	ITEM #	PRICE
Length 30cm (12in)	175490	\$8.10
Length 60cm (24in)	175493	\$8.60
Length 90cm (36in)	175495	\$8.60




## Fluid / Blood Delivery Set

DESCRIPTION	ITEM #	PRICE
Pressurize this set for high flow infusion. Includes six foot (180cm) tubing, drip chamber with filter, pressure vent, injection site and spin luer lock.	175496	\$12.35




## Arthroscopy Fluid Line Wound Flushing Line

DESCRIPTION	ITEM #	PRICE
High flow infusion set includes spike, ten feet (250cm) of tubing, clamp and spin luer lock. Fits basic peristaltic pump.	175499	\$6.60



## Plasmapheresis Catheter

DESCRIPTION	ITEM #	PRICE
12Ga x 13cm (5.25in) teflon catheter with 4 additional sideholes over 2.5cm of the distal end. Catheter-over-needle insertion method.	1212F	\$10.80



## Clear Cuff Pressure Infusor

High performance pressure infusor.

- Provides easy view of fluid bag and fluid level.
- Special built-in hook allows fluid bag to be replaced without removing pressure infusor from hanging device.
- Each unit includes the squeeze bulb assembly.


DESCRIPTION	ITEM #	PRICE
Sizes available		
1000ml	175500	\$46.85
3000ml	175510	\$116.40
5000ml		
Includes one 3000ml infusor and two extenders	175520	\$135.95
Inflation Squeeze Bulb Assembly		
Includes bulb and three-way stopcock.	175590	\$19.55



## Large Animal Trans-Tracheal Wash Kit / TTL

Catheter-over-needle insertion method. The semi-rigid flushing catheter has a stylet and luer adapter.

DESCRIPTION	ITEM #	PRICE
10Ga x 2.25in steel introduction catheter with 12Ga x 70cm (28in) flushing catheter for horses and foals	TW1228	\$23.55
14Ga x 2.75in teflon introduction catheter with 16Ga x 70cm (28 in) flushing catheter for large dogs and neonatal foals	TW1628	\$23.55



INTRODUCERS are catheter-over-needle

For small animal TTL and TTW, see page 23

## Broncho-Alveolar Lavage Catheter (BAL)

- Open end silicone heavy wall catheter for the lavage of selected segments of the lung fields using the nasotracheal approach.
- Distal balloon can be inflated to prevent retrograde fluid flow. 50cc balloon volume.
- Proximal two-way fitting allows simultaneous connection for both irrigation and aspiration.
- Reuseable.

DESCRIPTION	ITEM #	PRICE
10mmOD, 2.5mmID, 240cm (96in)	BAL240	\$72.75
10mmOD, 2.5mmID, 300cm (120in)	BAL300	\$72.75
Replacement Stopcock (10 pack)	2400	\$35.40

For small animal BALs, see page 23



## IV / Enteral Feeding Backpack and Surcingle

Use the IV BACKPACK for the mobile patient with a flat back. Allows free movement without the worry of lines being tangled. The pressure infusor allows for placement of a 1000ml fluid bag directly onto the patient's back. Variable infusion rates can be achieved by adjusting the fluid control dial. (Larger pressure infusors on page 31.)

The LINE MANAGEMENT SURCINGLE keeps lines from becoming tangled on a non-ambulatory patient.

DESCRIPTION	ITEM #	PRICE
IV Backpack	7000	\$121.55
IV Backpack Starter Kit Includes IV backpack, clear cuff pressure infusor, IV administration line and fluid control dial.	7001	\$181.30
Line Management Surcingle - 120cm (48in)	7002	\$41.20
Weanling Line Management Surcingle - 150cm (60in)	7005	\$45.00
IV Administration Line Includes drip chamber, pinch clamp, and rotating luer lock. Total length: 29in	7350	\$6.80
Fluid Control Dial	175480	\$8.90



IV Backpack Starter Kit




Line Management Surcingle

## Equine Tendon Sheath Flushing Catheter

- Provide localized infusion of medication and flushing directly into the affected area.
- Catheters fenestrated up to the 4cm or 8cm mark.
- 6Fr peel away introducer.

DESCRIPTION	ITEM #	PRICE
4cm of fenestrations	9120	\$28.00
8cm of fenestrations	9180	\$28.00



## Equine Joint Infusion Kit

Provide continuous localized infusion of antibiotics directly into the affected joint. A crush resistant 5Fr (16Ga) catheter inserted using a 6Fr (14ga) peel away introducer. A lightweight administration pump will hold up to 100ml and deliver 0.8ml/hr up to 2.5ml/hr.

ITEM # / PRICE	9100 / \$81.35
----------------	----------------



See MILA website for detailed insertion instructions written by Stephen B. Adams DVM, MS and Timothy B. Lescun, BVSc

## Cannula with Trocar

Introducer 10Ga x 5cm (2in) stainless steel cannula with trocar.

ITEM # / PRICE	1002 / \$12.00
----------------	----------------



## Accessories and Components

### Guidewire in Dispenser


One soft flexible 'J' end, one soft flexible straight end.

DESCRIPTION	ITEM #	PRICE
Sizes include:		
.018 x 45cm - Use with 4Fr & 5Fr catheters (MILA #1500, 1540, 1600, 1640, 1700, 1740, 1650, 1906, 1910, SA1606, SA1610)	GW1845	\$6.20
.032 x 45cm - Use with 16ga catheters (MILA #1606, 1606-2P, 1610-2P)	GW3245	\$6.20
.032 x 60cm - Use with 16ga catheter (MILA #1610, 1625-2P)	GW3260	\$6.20
.035 x 45cm - Use with 14ga and 7Fr catheters (MILA #1406, 1406-2P, 1410-2P, 1720)	GW3545	\$6.20
.035 x 60cm - Use with 14ga and 7Fr catheters (MILA #1410, 1620, 1630)	GW3560	\$6.20




### Vessel Dilator

DESCRIPTION	ITEM #	PRICE
4Fr x 5cm	DIL4	\$2.50
5Fr x 6cm	DIL5	\$2.50
6Fr x 6cm	DIL6	\$2.50
7Fr x 10cm	DIL7	\$2.50
8Fr x 10cm	DIL8	\$2.50
9Fr x 10cm	DIL9	\$2.50




### Double Christmas Tree Style Connector

DESCRIPTION	ITEM #	PRICE
Dual End 4-12mm (10 pack) Each individually packed and sterile	2300	\$32.75




### Double Christmas Tree Style Connector

DESCRIPTION	ITEM #	PRICE
4.0-6.0mm 7.5-11.5mm (10 pack) Each individually packed and sterile	2320	\$32.75




### Christmas Tree with Female Luer Connector

DESCRIPTION	ITEM #	PRICE
Adapter accepts 20-40Fr catheter (10 pack) Each individually packed and sterile	2350	\$32.75



### Christmas Tree with Male Luer Connector

DESCRIPTION	ITEM #	PRICE
Adapter accepts 20-40Fr catheter (10 pack) Each individually packed and sterile	2380	\$32.75



## Accessories and Components

### Stopcocks

DESCRIPTION	ITEM #	PRICE
Three-way with luer slip/10 pack Each individually packed and sterile	2400	\$35.40
Four-way with swivel luer lock/10pack Each individually packed and sterile	2410	\$35.40
Four-way with 18in extension	2415	\$7.20
Three-way with 18in extension	2418	\$7.20




### .2 Micron Filter

DESCRIPTION	ITEM #	PRICE
With male and female luer locks. (10 pack) Each individually packed and sterile	2420	\$47.00



### Female to Female Luer Connector

DESCRIPTION	ITEM #	PRICE
Connection of incompatible tubes/sets. (10 pack) Each individually packed and sterile	2440	\$18.55




### Male to Male Luer Connector

DESCRIPTION	ITEM #	PRICE
Connection of incompatible tubes/sets. (10 pack) Each individually packed and sterile	2445	\$18.55




### T Connector

DESCRIPTION	ITEM #	PRICE
Two female connectors with a male rotating luer lock. (10 pack) Each individually packed and sterile	2450	\$33.85




### Y Connector

DESCRIPTION	ITEM #	PRICE
Two female connectors with a male rotating luer lock. (10 pack) Each individually packed and sterile	2455	\$33.85



### W Connector


DESCRIPTION	ITEM #	PRICE
Three female connectors with a male rotating luer lock. (10 pack) Each individually packed and sterile	2460	\$33.85



## Accessories and Components

### Fluid Control Dial

DESCRIPTION	ITEM #	PRICE
Easily control infusion rates from 10ml/hour to 250ml/hour using gravity or 40ml/hr to 1000ml/hr using pressure fluid systems.	175480	\$8.90



### Suture Wings

DESCRIPTION	ITEM #	PRICE
10 per pack - individually packed and sterile		
4 French (Light Blue)	2530	\$29.85
5 French (White)	2550	\$17.35
7 French (Dark Blue)	2570	\$17.35



### Pinch Clamp

DESCRIPTION	ITEM #	PRICE
Small Pinch Clamp - 16Fr	Clamp -175420	\$6.70
Large Pinch Clamp - 20Fr	Clamp -175490	\$6.70

(15 per bag, non-sterile)



## Conversion Chart

FRENCH	GAUGE	INCH	MM
3	20	0.039	1
4	18	0.053	1.35
5	16	0.066	1.67
6	14	0.079	2
7	13	0.092	2.3
8	12	0.105	2.7
9	—	0.118	3
10	10	0.131	3.3
11	—	0.144	3.7
12	—	0.158	4
13	—	0.170	4.3
14	—	0.184	4.7
15	—	0.197	5
16	—	0.210	5.3
17	—	0.223	5.7
18	—	0.236	6
19	—	0.249	6.3
20	—	0.263	6.7
22	—	0.288	7.3
24	—	0.315	8
26	—	0.341	8.7
28	—	0.367	9.3
30	—	0.393	10
32	—	0.419	10.7
34	—	0.445	11.3

Length: 1 inch = 2.54 centimeters

## Ordering Information

1. Order direct or from a distributor.
2. Visa, MasterCard and American Express accepted.

### Mis-shipment

Notification must be given within 10 days of receipt of purchase for credit to be issued.

### Returns

1. Notification must be given before return.
2. Full credit will be allowed for the following conditions:
  - A. Merchandise shipped in error against purchase order.
  - B. Defective merchandise.
3. When unopened merchandise is authorized to be returned for any reason other than above, credit will be:
  - 80% - less than one year old
  - 50% - one to two years old
  - No credit - more than two years old

Prices Subject to Change Without Notice

ITEM #	PAGE										
1006	9	3240	11	175497	30	EMAC700L	28	J636	18	PI-1740	7
1211	9	3360	11	175498	30	EMAC800	28	J836	18	PI-1906	7
1212	9	3460	11	175499	31	ETC1033	25	JP4	22	PI-1910	7
1212F	31	3560	11	175500	31	ETC1055	25	JP7	22	PI-506	7
1402	9	3565	11	175510	31	ETC1133	25	JP10	22	PI-510	7
1403	9	4203	8	175520	31	ETC1155	25	JP13	22	RAFIX16	18
1406	6	4204	8	175590	31	ETC1233	25	LAN190	28	RAFIX20	18
1406-2P	6	4205	8	175720	30	ETC1255	25	LL1945	10	RPEG16	18
1410	6	4205.5	8	175800	21	ETC20	25	LL1990	10	RPEG20	18
1410-2P	6	4206	8	175810	21	ETC30	25	LL2045	10	SA1606	6
1411	9	4206.5	8	175820	21	ETC40	25	M555	17	SA1610	6
1500	6	6612	27	175830	21	ETC50	25	MC40	12	SF10	26
1540	6	6612L	27	8V0002	30	ETC755	25	MC43	12	SYR10	26
1600	6	6613	27	A1410	6	ETC855	25	MC44	12	T1803	19
1603	9	6613L	27	A1610	6	ETC933	25	NG1036	16	T1806	19
1606	6	6615	27	A1620	6	ETC955	25	NG1036W	16	T2025	19
1606-2P	6	7000	32	API-1410	7	ETN10	25	NG1243	16	T2035	19
1610	6	7001	32	API-1610	7	ETN15	25	NG1243W	16	TW1228	32
1610-2P	6	7002	32	BAL240	32	ETN20	25	NG1450	16	TW1610	23
1611	9	7005	32	BAL300	32	ETN25	25	NG18100	16	TW1616	23
1620	6	7100	26	BL8	22	ETN30	25	NG2496	16	TW1617	23
1625	6	7350	32	BL11	22	ETN35	25	NG522S	16	TW1628	23/32
1625-2P	6	7360	32	BMAN	19	ETN40	25	NG622	16	UC1000	15
1630	6	8050	13	CF34110	39	FCT REV	26	NG622W	16	UC1000K	15
1640	6	8055	13	CFS36	19	FCT12	26	NG636	16	UC1018	15
1650	6	8098	13	Clamps	36	FCT15	26	NG843	16	UC1036	15
1660	6	8099	13	CT1215	20	FCT2	26	NG843W	16	UC1212	15
1700		8080	13	CT1410	20	FCT30	26	NG855	16	UC1221	15
1720	6	8090	13	CT1415	20	FCT4	26	NG1043	16	UC1621	15
1740	6	8095	13	CT2015	20	FIX16	18	NG1055	16	UC149	15
1760	6	8515	12	CT2022	20	FIX20	18	NL325	27	UC2416	15
1806	19	8518	12	CT3015	20	G100	22	NL525	27	UC2459	15
1808	19	8520	12	CT3022	20	G150	22	NL555	27	UC306	15
1810	19	8540	12	Cubby	17	G200	22	NL525FS	27	UC306S	15
1811	9	8570	12	CVP35	11	G400	22	OCN16	29	UCC308	15
1820	6	8572	12	D1090	22	GDV1	23	OCRL	29	UCO308	15
1906	6	8575	12	D890	22	GDV2	23	OCRN	29	UC310	14
1910	6	8585	12	DIFF2	24	GDV3	23	OCRN-2	29	UC310S	14
2008	19	8590	12	DIFF4	24	GW1845	34	OCS	29	UC506	14
2011	9	8599	12	DIFF6	24	GW3245	34	OCS-2	29	UC506S	14
2206	19	9100	33	DIFF7	24	GW3260	34	OCSA	29	UC510	14
2217	9	9120	33	DIFF9	24	GW3545	34	PEG16	18	UC510S	14
2300	34	9180	33	DIL4	34	GW3560	34	PEG20	18	UC624	15
2320	34	9203	24	DIL5	34	H2900	21	PGS16	18	UC636	15
2350	34	9206	24	DIL6	34	IC1801	10	PGS20	18	UC812	15
2380	34	9500	28	DIL7	34	IC2001	10	PI-1410	7	UC824	15
2400	35	9510	28	DIL8	34	IC2201	10	PI-1500	7	UC836	15
2410	35	9520	28	DIL9	34	IC2475	10	PI-1540	7	UFC3480	29
2415	35	175420	21	DK1090	22	ID010	27	PI-1600	7	UFC4080	29
2418	35	175425	21	DK890	22	ID06	27	PI-1606	7	UFYA	29
2420	35	175440	21	DVD100	13	ID08	27	PI-1606-2P	7	V2030	11
2440	35	175450	21	E1413	17	ID10	27	PI-1610	7	WIC1801	10
2445	35	175460	21	E1430	17	ID14	27	PI-1610-2P	7	WIC2001	10
2450	36	175470	21	E1450	17	ID21	27	PI-1624	7	WIC2201	10
2455	36	175480	35	E1830	17	ID42	27	PI-1625	7	WIC2475	10
2460	36	175488	30	EDC1624	23	ID50	27	PI-1625-2P	7		
2530	36	175490	31	EDC1648	23	ID-T	27	PI-1640	7		
2550	36	175493	31	EDC190	28	ID-W	27	PI-1648	7		
2570	36	175495	31	EDC400	28	INS200	28	PI-1650	7		
		175496	30	EMAC700	28	J1036	18	PI-1700	7		

**AUSTRALIA**

CPP Cardio Pulmonary Products  
Leichhardt  
email: cpp112@bigpond.net.au  
Phone: 61 02 9569 6810

Vetquip  
NSW  
email: vetquip@bigpond.com  
Phone: 61 02 9889 6666

**BELGIUM & THE NETHERLANDS**

Instrulife NV  
Oostkamp  
email: info@instrulife.be  
www.instrulife.be  
Phone: +32 50 311876

**BRAZIL**

Digicare Animal Health  
Rio de Janeiro  
email: vendas@digicarebiomedical.com  
www.digi-vet.com.br  
Phone: +21-3373-1378

**CANADA**

Central Sales Limited  
Brampton, ON  
800-387-2522

Pacific Veterinary Sales Ltd  
Abbotsford, BC  
Phone: 604-850-3555  
Fax: 604-850-3530

**CHINA**

ANSCITECH Limited  
Wuhan, Hubei  
email: info@anscitech.com  
www.anscitech.com  
Phone: +86 27 50705661

**COLOMBIA**

VitaeQuinos  
Medellin, Colombia  
www.vitaequinos.com  
Phone: +574 3116799

**DENMARK**

E-vet  
Haderslev  
email: salg@evet.dk  
www.evet.com  
Phone: +45 7453 0848

**FINLAND**

Vetman  
Halikko  
Phone: 358 2 7221 222

**FRANCE**

ALCYON  
Paris  
email: contact@alcyon.com  
Phone: 3301 53903939

Coveto  
Montaigu  
email: materiel@coveto.fr  
Phone: 33 02 51488088

Genia  
St. Hilaire de Chaleons  
email: H.Chartier@genia.fr  
Phone: 33 02 4033 3537

**GERMANY-AUSTRIA-SWITZERLAND**

BioRanch  
Ostbevern  
email: andre.schaefer@bioranch.de  
Phone: 49 2532 959 385

**HUNGARY**

Hungarovet Allatkorhaz Kft.  
Budapest  
www.hungarovetproducts.com  
Phone: 36 1 362 2699

**ISRAEL**

Dyn Medical Equipment, Ltd.  
www.dyn.co.il  
Phone: +972 4 6175331

VetMarket  
Hod Hasharon  
email: vetmarkt@vetmarket.co.il  
Phone: 972 9-7414145

**ITALY**

Alcyon Italia Spa  
Marene - CN  
email: marketing@alcyonitalia.com  
Web: www.alcyonitalia.com  
Phone: +39 0172 743111

**JAPAN**

AVS Co., Ltd  
Tokyo  
www.av.s.cojp/en/  
Phone: 81 3 3340 2014

**MEXICO**

Difarvet S.A. de C.V.  
email: abenavides@difarvet.com  
www.difarvet.com  
Phone: +52 33 3836-0070  
52 33 800-024-0928

**NETHERLANDS**

Instrulife NV  
Oostkamp  
email: info@instrulife.be  
Phone: +32 50 311876

**NEW ZEALAND**

MEDSOURCE  
Ashburton  
email: info@medsource.co.nz  
Phone: 64 33089718

SHOOF International Ltd.  
Cambridge  
email: sales@shoof.co.nz  
Phone: 64 7 827 3902

SVS Veterinary Supplies Ltd.  
Hamilton  
email: sales@svs.co.nz  
Phone: 64 39 02 642 8299

**NORWAY**

Jan F. Andersen A/s  
Jevnaker  
email: sales@jfa-as.no  
Phone: 47 61 31 4949

**PORTUGAL**

Vetin  
Quinta do Anjo  
www.vetin@sapo.pt  
Phone: 21 213 7660

**SPAIN**

Ekufarm S.L.  
Barcelona  
Phone: 34 938 459400

LaBouvet  
Torrejon de Ardoz  
email: labouvet@labouvet.com  
Phone: 34 91 726 4229

**SWEDEN**

SweVet Piab AB  
Sjoberg  
email: info@swevet.se  
Phone: 46 416 258 27

Vetgross  
Blentarp  
email: info@pharmadistribution.se  
Phone: 46 416 246 60

**UNITED KINGDOM**

Direct Medical Supplies, Ltd.  
Hampshire www.dmsveterinary.com  
email: sales@directmedical.uk.com  
Local Phone: 7041 527 007  
Local Fax: 7041 527 008  
Phone: 44 1420 520 869

## **Butler Animal Health Supply Company**

---

Dublin, OH  
800-551-3861

## **DVM Resources**

---

Grapevine, TX  
877-828-1026

## **Lextron Supply**

---

Lago Vista, TX  
888-787-4483

## **Midwest Veterinary Supply**

---

Burnsville, MN  
800-352-2831

## **MWI Veterinary Supply Co.**

---

Nampa, ID  
800-824-3703

## **Penn Veterinary Supply, Inc.**

---

Lancaster, PA  
800-233-0210

## **Webster Veterinary Supply**

---

Sterling, MA  
800-225-7911



**MILA**   
**INTERNATIONAL, INC.**  
*Medical Instrumentation for Animals*

Your MILA Distributor:

Order direct

MILA International, Inc.  
12 Price Avenue  
Erlanger, KY 41018 USA

Phone +1-859-957-1722  
FAX +1-859-957-1727  
Toll free 888-MILAIN (645-2468)  
Web [www.milaint.com](http://www.milaint.com)